

# TINYSCOPE

## PRODUCT MANUAL

TinyScope

Thank you for purchasing and using

### Product introduction

#### 1. Packing list

TinyScope microscope, spare sticker, lanyard, lens wipe, storage box, instruction, glass slide.

#### 2. Specification



### Manual

1. Download the APP:  
IOS: Search "TinyScope" in the Store.  
Android: Search "TinyScope" in Google Play.  
(ps: System camera APP can work if the flash can keep lighting)

2. Open the storage box by pulling the corner of the buckle, and take out the TinyScope.




3. Put the TinyScope facing the table, and open the flash in the APP.



4. Take off the phone shell, and aim the camera of the mobile phone to the lens of TinyScope, while the flash of the mobile phone facing the body of TinyScope. Move the mobile phone to make the field of view in the center of the screen.



5. Move the mobile phone towards the TinyScope until touching. Push the mobile to make TinyScope stick to the mobile.



6. Move TinyScope close to the object.



(ps: TinyScope fits most mobile phones)

### Replace sticker

1. Take off the old sticker, and get a new spare sticker in the storage box.



2. Get one piece of spare sticker, and tear off the release paper of one side. Please do not touch the surface of the sticker.



3. Stick the sticker onto the TinyScope, and press.



4. Tear off the release paper of the other side.



### Product parameters

Product brand	TinyScope
name	TinyScope microscope
Model	TS-V1
Material	PC
Size	42 mm×18 mm×5 mm
weight	5 g

### Notice

1. TinyScope is not waterproof. TinyScope is battery free, uses mobile flash for illumination, and is stick-and-play. If there is more than one mobile camera, please make the lens in front of the main camera.

2. Do not touch the sticker of TinyScope, and keep it clean.

3. Please keep the lens away from scratch and pollution. Please use lens wipe to clean. Please keep TinyScope away from wet. We suggest put TinyScope into the storage after playing.

4. The size of TinyScope is small, please keep it unavailable for children.

### Exemption Instruction

Thank you for purchasing and playing TinyScope.

Out of warranty: damaged by wet, dirt, scratching.

If the cover glass of the mobile phone camera drops, it's the mobile phone manufacturer's blame, and our company has no responsibility.

If there is any differences between instruction and product, please take the product as the standard.

If you have any question or objection on this instruction, please write to us in seven days after purchasing. Or it is considered that you understand and agree to the entire instruction. Our company keep the right to explain and modify the instruction within the law.

### Exchange policy

1. Please keep the original appearance of product when refund or exchange (except quality fault of product during warranty period).

2. Original appearance of product includes TinyScope microscope, spare sticker, sling, lens cloth, storage box, and instruction. Appearance should not affect later sale.

3. Sellers pay the freight fee of refund/exchange if it is caused by quality fault of product.

4. Customers pay the freight fee of refund/exchange if it is not caused by quality fault of product.

5. Call 86-027-65527110 if you want refund/exchange.

Manufacturer: Convergence (Wuhan) Technology Ltd., China

Tele hone number: 86-027-65527110

Address: 2505 Huigushikong building, Luoyu Road, Hongshan District, Wuhan, Hubei, China

Website: www.tinyscope.com